# i-BONDING® LC<sup>N</sup>

EN	Light curing adhesive
LT	Šviesoje kietėjantis surišėjas
DE	Lichthärtender Klebstoff
BG	Леко вътвърдяващо лепило
CS	Světlem tuhnoucí lepidlo
DA	Let hærdende klæbermiddel
ET	Kergelt kövenev liim
FR	Adhésif photopolymérisable
EL	Κόλλα ελαφριάς σκλήρυνσης
HU	Fényre kötő ragasztó
IT	Adesivo fotopolimerizzabile
LV	Gaismas cietēšanas līme
NO	Lettherdende lim
PL	Klej światłoutwardzalny
PT	Adesivo fotopolimerizável
RO	Adeziv de întărire usoară
SK	Svetlo tuhnúce lepidlo
SL	Lepilo za utrjevanje svetlobe
ES	Adhesivo fotopolimerizable
SV	Lätthärdande lim
TR	Işıkla sertleşen yapıştırıcı
RU	Светоотверждаемый клей
UA	Світлотвердючий клей

EN

## INSTRUCTION FOR USE

### DESCRIPTION

i-BONDING LC<sup>N</sup> is a 5<sup>th</sup> generation light curing, bioinert, single component one-step adhesive, suitable for use with wet bonding technique. Etching of the surfaces is required.

### COMPOSITION

Methacrylate mixture 55-95%, acetone 10-30%, conditioner 1-10%, coinitiator <5%, accelerator <1%, inhibitor <1%, photoinitiator <1%. i-BONDING LC<sup>N</sup> does not contain medicinal substance, including human blood or plasma derivate; tissues or cells, or their derivates, of human origin; tissues or cells of animal origin, or their derivates, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

### PERFORMANCE CHARACTERISTICS

light curing	20 s
shear bond strength (enamel)	21.8±2.4 MPa

### INTENDED PURPOSE AND CLINICAL BENEFITS

i-BONDING LC<sup>N</sup> helps to restore/improve aesthetic appearance of restorable tooth; helps to restore/maintain dental function of restorable tooth; helps to protect biological structures of restorable tooth and adjacent tissues.

### CLINICAL INDICATIONS

- Used to bond dentin and enamel to polymer-based restorative materials, compomers, ceramic surfaces to improve retention of a restoration.

### CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

### RESTRICTIONS TO COMBINATIONS

i-BONDING LC<sup>N</sup> should not be used with products containing eugenol as eugenol may disturb polymerization process.

### UNDESIRABLE SIDE EFFECTS

In susceptible individuals, i-BONDING LC<sup>N</sup> may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract).

### RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

### PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

### INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth.

### INTENDED USER

i-BONDING LC<sup>N</sup> is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental adhesive. There is no need for specific training.

### STERILITY

i-BONDING LC<sup>N</sup> is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

### USE ENVIRONMENT

i-BONDING LC<sup>N</sup> is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of adhesive is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

### CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as micro applicators, cup for application, are supplied with the device.

### INSTRUCTION FOR USE

1. Prepare cavity as always. Cleaning with oil-free prophylaxis paste, such as i-FASTE is recommended. For deep cavities use calcium hydroxide liner or glass ionomer based cement, such as i-BAS.
2. Apply layer of etch, such as i-GEL<sup>N</sup> on all surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin. After drying, the surface must be chalky and do not contaminated prior to applying the bonding resin. If the surface becomes contaminated with saliva, rinse and dry again and re-etch if necessary.
3. Shake bottle before each use.
4. Delicate push on bottle and dispense small amount of adhesive into cup. Do not force bottle.
5. Using a fully saturated micro applicator apply 2-3 layers of adhesive onto etched surface.
6. Dry with air for 10-20 seconds to dry all solvent from surface. If surface is not glossy, apply additional layers of adhesive and dry.
7. Curing with light for 20 seconds. Use LED polymerization lamp with light intensity 1200mW/cm<sup>2</sup>. Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use.
8. Apply the layer of light curing restorative material, such as i-LIGHT<sup>N</sup>, i-XCITE LC<sup>N</sup> or i-FLOW<sup>N</sup> into cavity and finish restoration.

### WARNINGS

After the desired amount of adhesive extruded, immediately close the bottle. i-BONDING LC<sup>N</sup> is highly flammable liquid. Do not use i-BONDING LC<sup>N</sup> for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. i-BONDING LC<sup>N</sup> does not emit radiation and does not cause any electromagnetic interferences.

### PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing.

Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

### SHELF-LIFE

Shelf-life of i-BONDING LC<sup>N</sup> is 3 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

### STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

### DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

### VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

i-BONDING LC<sup>N</sup> is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

### MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

### VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

### PACKAGING

REF IBLTP 5ml bottle, 5pcs micro applicators, cup

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